



ACCUTRON[®] CT-D VISION GANTRY VERSION FOR SIEMENS HEALTHINEERS

MEDTRON presents its CT scanner mounted version of the dual head injector Accutron[®] CT-D Vision.

It offers the best of MEDTRON's dual head injection technology & benefits to SOMATOM go CT scanners with a physical integration onto the scanner.

**MEDTRON AG -
ESSENTIAL FOR CONTRAST**



ULTIMATE FLEXIBILITY IN CT IMAGING

The new Accutron® CT-D Vision Gantry Version provides the ultimate flexibility in contrast-enhanced CT imaging when using a power injector.

MEDTRON AG in cooperation with Siemens Healthineers has designed a natural extension of the SOMATOM go CT scanners, providing maximum injector uptime with complete Accutron® CT-D Vision functionality. The Accutron® CT-D Vision Gantry Version (CT861-2), a dual head injector directly mounted on the gantry by means of support arm by Siemens Healthineers ensures free range of operations in the examination room and a closer-to-patient positioning on either side of the bore to deliver the full spectrum of contrast enhanced procedures.


FULL INTEGRATION

WITH SIEMENS HEALTHINEERS
CT-SCANNERS

Combining Accutron® CT-D Vision and SOMATOM go CT scanner with MEDTRON's CANopen Class IV interface allows synchronization to achieve an optimal workflow. Thus Accutron® CT-D Vision Gantry Version provides full integration with the CT scanners.

The CanOpen Class IV interface for CT scanners from Siemens Healthineers (IF 864) is available as an option.

Accutron® CT-D Vision Gantry Version is delivered with the following accessories:



- 1 Contrast / Saline containers suspension hooks
- 2 Permanent power supply is ensured independently through the gantry support arm to the SOMATOM go CT scanner. (The gantry support arm is supplied by Siemens Healthineers.)
- 3 Hand switch to launch the injection
- 4 Side arm with cup to vent the tube system
- 5 Wireless touch screen remote control

The Accutron® CT-D Vision Gantry Version is compatible with :

SOMATOM go.Now, SOMATOM go.Up,
SOMATOM go.All, SOMATOM go.Top.

Installation note: The gantry support arm by Siemens Healthineers is required to install the injector.

Clinical Gains



- The optimal integration between Accutron® CT-D Vision and the gantry support arm by Siemens Healthineers provides an easier to clean environment to maintain and enhance proper hygiene protocols.

Operating Gains



- This integration provides close patient positioning of a power injector whether it is positioned near the patient's head, or opposite of the CT bore with proximal positioning to the lower extremities. This ensures an efficient procedure of the imaging examination.
- Accutron® CT-D Vision Gantry Version adds to ones flexibility. Due to constant power supply and low maintenance effort, it provides maximum uptime with complete Accutron® CT-D Vision functionality.
- Mounted to the CT scanner, the injector offers free range of operations in the examination room. The permanent power supply is provided by a cable in the gantry support arm to the SOMATOM go CT scanner.
- For an optimal workflow the injector, equipped with an optional CANopen Class IV interface, allows synchronization to achieve full integration with Siemens Healthineers CT Scanners.

Financial Gains



- By providing a large scope of applications, increasing patient workflow and simplifying maintenance, Accutron® CT-D Vision Gantry Version contributes to a better return on investment.
- Even as an accessory to entry level CT scanners like SOMATOM go, Accutron® CT-D Vision offers all the functionalities one can expect from an advanced dual head contrast media injector.

Please contact your MEDTRON representative or refer to your according sales contact at Siemens Healthineers in your country, if you wish to obtain

- the gantry support arm only
- the gantry support arm with injector
- a CT scanner with gantry support arm and injector as full solution

MEDTRON AG | Hauptstrasse 255 | 66128 Saarbruecken, Germany
Tel.: +49 (0)681-97017-0 | Fax: +49 (0)681-97017-20 | info@medtron.com

Visit us online at www.medtron.com



Contact:

Team DACH:
Germany, Austria, Switzerland
Tel.: +49 (0)681-97017-72
Fax: +49 (0)681-97017-60
sales.dach@medtron.com

Team International 1:
W/S-EMEA, LATAM, Africa, APAC
Tel.: +49 (0)681-97017-26
Fax: +49 (0)681-97017-20
sales.int1@medtron.com

Team International 2:
E-Europe, CIS
Tel.: +49 (0)681-97017-63
Fax: +49 (0)681-97017-20
sales.int2@medtron.com