



MED (TRON[®]) AG

Efficient
&
easy to use

The flexible solution!

Consumables - newly approved!
For multiple use up to 8 hours.

Advantages and new aspects of our multiple use consumables



YOUR BENEFITS

AS COMPARED TO DISPOSABLES:


- ✓ Economically efficient by time and cost savings!
- ✓ Can be used up to 8 hours!
- ✓ Overview and control by special stickers!
- ✓ Individual system solution!
- ✓ Simple handling!
- ✓ Environmentally compatible!
- ✓ Compact and well-structured design!

IT'S THIS SIMPLE:

You can save time and material with our customisable modular system comprising of an 8h-syringe, 8h-tube system and patient line (with green check valve).

Look forward to up to 8 hours without syringe and tube system changes. During this period, the products can be used on different patients. You only need to swap the patient line after every patient.*

We do our best for a simple handling of our products. Our sophisticated system is user-friendly and will significantly simplify your workflow during the examination!



Data and facts

Syringe

- Article number 31 4626
ELS 200 ml syringe MU

Selection of all approved tube systems and patient lines for 8h multiple use:


Tube system with yellow check valve

- Article number 31 4081
EK 223 tube system with drip chamber
- Article number 31 4099
MRS 222 MR-tube system with inline drip chamber
- Article number 31 4100
MRS 222 MR-tube system
- Article number 31 4108
Double head tube system with inline drip chamber and valve
- Article number 31 4119
MRS 222 XS tube system with replaceable spike and additional valve

Patient line with at least one green check valve

- Article number 31 8020
ES 224/25 Patient line with valve
- Article number 31 8026
Patient line 150 cm with terminal check valves
- Article number 31 8036
Patient line 150 cm with 2 valves
- Article number 31 8152
Patient line 150 cm with valve
- Article number 31 8181
Coiled line 180 cm with valve
- Article number 31 8198
Patient line 200 cm with valve
- Article number 31 8214
Coiled line 150 cm with valve
- Article number 31 8217
Coiled line 150 cm with 2 valves
- Article number 31 8251
Patient line 250 cm with valve
- Article number 31 8301
Patient line 300 cm with valve

The syringe is approved in combination with any of the tube systems for multiple use within a duration of up to 8 hours!



Declaration of Compatibility

Hygienic safety of the consumables for 8h-application has been confirmed by an independent test lab.

MEDTRON[®] AG

MEDTRON AG · Hauptstr. 255 · 66128 Saarbrücken
Phone +49 (0) 681 97017-25 · Fax +49 (0) 681 97017-20
Email info@medtron.com · www.medtron.com



* Please refer to our manual, the general hygiene principles and the special hygiene requirements of your company in order to prevent the product from being contaminated from the outside. Also, please be sure to adhere to the indicated maximum utilisation periods, including those of the contrast medium and the sodium chloride solution (see package insert)